



Key Features

- Self-adjusting flexible forehead support allows for a wide range of angle positions
- Position of forehead support can be adjusted vertically to three different heights
- Easy one hand adjustment is perfect for patients with limited dexterity
- Quick release clips
- 360° rotation of the elbow
- Perfectly contoured soft, lightweight silicone cushion
- Cushion inflates to facilitate sealing therefore requiring less headgear tightening
- 100% latex free

SYLENT Silicone NIV Nasal Mask

The Sylent NE Disposable Mask

Quick, comfortable, quality

The Sylent NE disposable nasal mask is designed for providing non-invasive ventilation or pressure support therapy. The completely flexible silicone forehead support offers infinite adjustability and quick release clips allow the mask to be removed and worn again without adjustment.

SIZES

Small

Medium

MATERIAL

Cushion	Silicone
Caps	Silicone
Shell	Polycarbonate
Elbow	Polycarbonate
Rings	Polycarbonate
Forehead Support	Polycarbonate
Buckles	Polypropylene
Headgear	Nylon, Neoprene, Velcro Hook, Loop



ORDERING INFORMATION

RespCare nr.	Model	Size	Order nr.
SYL271	Sylent NE Disposable Mask	Small	417 005 000
SYL273	Sylent NE Disposable Mask	Medium	417 005 001

Manufactured for

RESPCARE



US: Rx only

Distributed by

Ambu A/S

Baltorpbakken 13

DK-2750 Ballerup

Danmark

Tlf.: +45 72 25 20 00

Fax: +45 72 25 20 53

www.ambu.com